TA Instruments

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Q Series™ is a trademark of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

Integrity™ is a trademark of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

Modulated DSC® and MDSC® are registered trademarks of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720

Tzero™ is a trademark of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

μTA® is a registered trademark of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

Smart Swap™ is a trademark of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

Hi-Res™ is a trademark of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

Mobius Drive™ is a trademark of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

Modulated TGA® and MTGA® are registered trademarks of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

Platinum™ is a trademark of TA Instruments Waters—LLC, 109 Lukens Drive, New Castle, DE 19720.

Patente von TA Instruments

Apparatus for Controlling the Power Supplied to a Load describes proprietary technology patented by TA Instruments Waters—LLC (CA Patent No. 1,133,988).

Method and Apparatus for Calorimetric Differential Thermal Analysis describes proprietary technology patented by TA Instruments Waters—LLC (CA Patent No. 1,159,278).

Autosampler describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 4,816,730).

Parallel Plate Dielectric Analyzer describes proprietary technology patented by TA Instruments Waters—LLC (EP Patent No. 0347125).

Method and Apparatus for Modulated Differential Analysis (MDSC®) describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent Nos. 5,224,775; 5,248,199; 5,346,306; 5,439,291. Additional Patent Nos. CA 2,089,225; JP 2,966,691 and 3,299,575.

Heat Flux Differential Scanning Calorimeter Sensor (Tzero™) describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent Nos. 6,431,747 and 6,561,692).

Power Compensation Differential Scanning Calorimeter describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,428,203).

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Method and Apparatus of Modulated-Temperature Thermogravimetry (MTGA™) describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent Nos. 6,336,741 and 6,113,261).

Infrared Heated Differential Thermal Analyzer describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 5,509,733. Additional Patent Nos. CA 2,134,432 and BE, DE, EP, FR, and GB 0660110).

Mechanical Cooling System describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 5,484,204. Additional Patent Nos. BE, CH, DE, EP, FR, GB, LI, and NL 0703448).

Method and Apparatus for Performing Localized Thermal and Sub-Surface Imaging by Scanning Thermal Microscopy describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent Nos. 6,095,679 and 6,491,425. Additional Patent No. JP 3,229,329 and CA 2,225,315).

Modulated-Temperature Thermomechanical Analysis describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,007,240).

Method and Apparatus for High Resolution Analysis of the Composition of a Material describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 5,368,391 and 5,165,792. Additional Patent Nos. CA 2,051,578 and DE, EP, FR, GB, IT 0494492).

Method and Apparatus for Thermal Conductivity Measurements describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 5,335,993).

Dynamic and Thermal Mechanical Analyzer Having an Optical Encoder with Diffraction Grating and a Linear Permanent Magnet Motor describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent Nos. 5,710,426 and 5,915,283. Additional Patent No. EP 0793088).

Heater Control Evaporation of Cryogenic Fluids for Cooling Thermal Analysis Instruments describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 5,973,299).

Differential Scanning Calorimeter describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 5,842,788 and 6,641,300).

Method and Apparatus for Performing Chemical Analysis Using Imaging by Scanning Thermal Microscopy describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 6,405,137).

Thermogravimetric Apparatus describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 5,321,719).

Method and Apparatus for High Spatial Resolution Spectroscopic Microscopy describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 6,260,997).

Method and Apparatus for Localized Mechano-Thermal Analysis Combined with Scanning Probe Microscopy describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 6,200,022 and GB Patent No. 2,332,949).

Differential Scanning Calorimeter describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 6,488,406).

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Thermal Analysis Assembly with Distributed Resistance and Integral Flange for Mounting Various Cooling Devices describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 6,523,998).

Liquid Nitrogen Cooling System describes proprietary technology patented by TA Instruments Waters—LLC (U.S. Patent No. 6,578,367).

Rotary Rheometer (Peltier Concentric Cylinder) describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,588,254).

Method and Apparatus for Spatially Resolved Modulated Differential Analysis describes proprietary technology patented by TA Instruments—Waters LLC (EP Patent No. 0572164).

Thermopile Differential Thermal Analysis Sensor describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 5,288,147).

Gripper Device describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,652,015).

Method and Apparatus for Positional Calibration of a Thermal Analysis Instrument describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,760,679 and 6,823,278).

Sample Tray for Automatic Sampler describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,644,136).

System and Method for Calibrating Contact Thermal Resistances in Differential Scanning Calorimeters describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,648,504).

Apparatus and Method for the Simultaneous Measurement of Rheological and Thermal Characteristics of Materials and Measuring Cell describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 5,520,042).

Automatic Sampling Device (Q1000/100 Autosampler) describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,846,455).

Platen for Automatic Sampler (Q1000/100 Autosampler) describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,859,271).

Differential Scanning Calorimeter Accounting for Heat Leakage describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,843,595).

Devices, Systems and Methods for Sensing Temperature of a Drag Cup in a Rheometer Motor describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,798,099).

Thermogravimetric Analyzer Autosampler Sealed Sample Pan (Q5000) describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,840,668).

Wide Range Dynamic Rheometer describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,931,915).

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System and Method for Automatic Identification of a Detachable Component of an Instrument describes proprietary technology patented by TA Instruments - Waters LLC (U.S. Patent No. 6,952,950).

AR-G2 Smart Swap Geometries describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 6,952,950).

Rotary Rheometer Magnetic Bearing describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 7,017,393)

Differential Scanning Calorimeter Accounting for Heat Leakage describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 7,025,497)

Humidity Controlled Chamber for a Thermogravimetric Instrument describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 7,048,435)

System and Method for Measurement of Small-Angle or Small Displacement describes proprietary technology patented by TA Instruments—Waters LLC (U.S. Patent No. 7,075,317)

Polymer Melt and Elastomer Extension Fixture describes proprietary technology patented by TA Instruments - Waters LLC (U.S. Patent No. 7,096,728)

Absolute Temperature Measuring Apparatus and Method describes proprietary technology patented by TA Instruments - Waters LLC (U.S. Patent No. 6,994,467)

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Windows® NT, 2000, XP, 98, 98SE, Me, Microsoft Excel® and Microsoft Word 97® are registered trademarks of the Microsoft Corporation.

Adobe® Acrobat® Reader® are registered trademarks of Adobe Systems Incorporated.

Oracle® and Oracle9i™ are trademarks or registered trademarks of Oracle Corporation.

TrueMetrix™ and Scanning Tip Technology™ are registered trademarks of ThermoMicroscopes, Inc.

CHROMEL® and ALUMEL® are registered trademarks of Concept Alloys.

Teflon® is a registered trademark of E. I. du Pont de Nemours and Company.

Loctite® is a registered trademark of the Loctite Corporation.

Swagelok® is a registered trademark of the Swagelok Company.

Inconel® is a registered trademark of Inco Alloys/Special Metals.

X-acto® is a registered trademark of Hunt Corporation.

Novacure® and Omnicure™ are trademarks of EXFO Photonic Solutions, Inc.

Kalrez® is a registered trademark of DuPont Dow Elastomers.

TA Instruments Q Series modules contain proprietary embedded operating system software copyrighted by Mentor Graphics.

SILICON SOFTWARE

©1989-97 Mentor Graphics Corporation, Microtec Division. All rights reserved. Unpublished-rights reserved under the copyright laws of the United States.

RESTRICTED RIGHTS LEGEND

Use duplication, or disclosure by the U.S. Government or a U.S. Government subcontractor is subject to the restrictions set forth in the license agreement provided with the Software pursuant to DFARS 227.7202-3(a) or as set forth in subparagraph (c) (1) and (2) of the Commercial Computer Software-Restricted Rights clause at FAR 52.227-19 as applicable.

MENTOR GRAPHICS CORPORATION, MICROTEC DIVISION,
880 RIDDER PARK DRIVE, SAN JOSE, CA 95131-2440

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Email: INFO@tainstruments.com

HELPLINE—U.S.A.

For assistance with thermal analysis applications,
please call the Thermal Analysis Help Desk
at 1-302-427-4070.

SERVICE—U.S.A.

For instrument service and
repairs, please call 1-302-427-4050.

AUSTRALIA

AUSTRALIA

TA Instruments

C/O Waters Australia Pty. Ltd.

Unit 3, 38-46 South Street

Rydalmere NSW2116

Steve Shamis

Phone: 61-3-9553-0813

Fax: 61-3-9553-0813

BELGIUM / LUXEMBOURG

TA Instruments - Division of Waters NV/SA

ConneXion Business Park

Brusselsesteenweg 500

B-1731 Zellik

Belgium

Phone: 32/2 706 00 80

Fax: 32/2 706 00 81

E-Mail: Belgium@tainstruments.com

CANADA

c/o TA Instruments – Waters LLC

270 Inglehart Road

RR #1 Grimsby, Ontario

L3M 4E7 Canada

Phone: 905-309-5387

Fax: 905-309-5492

E-mail: shunt@tainstruments.com

CHINA

TA Instruments—Waters LLC

Hotline: 800-820-3812

Shanghai:

Waters International Trading

(Shanghai) Limited

16 F, No. 82 Building,

No. 1198, Qin Zhou Rd.(N)

Cao He Jing Developing Park

Shanghai, P.R. CHINA 200233

Phone: 8621 54263960

Fax: 8621 64956366

Email: info@tainstruments.com.cn

Website: <http://www.tainstruments.com.cn>

Beijing:

22F, No.3 Building, No. 98, West

Ba Li Zhuang, Chaoyang District,

Beijing, 100025

Phone: 8610 85868899

Fax: 8610 85867099

E-mail: info@tainstruments.com.cn

Website: <http://www.tainstruments.com.cn>

Guangzhou:

Rm. 406-407, Commercial Building, China Hotel

Liuhua Road

Guangzhou 510015

Phone: 8620 86266678

Fax: 8620 86686217

Email: info@tainstruments.com.cn

Website: <http://www.tainstruments.com.cn>

Hong Kong:

Rm. 901, 102 Austin Road

Kowloon

Hong Kong

Phone: 852 29641800

Fax: 852 25496802

Email: info@tainstruments.com.cn

Website: <http://www.tainstruments.com.cn>

FINLAND

Waters Finland Oy

Kuparitie 1

Fin-00440 Helsinki

Finland

Phone: 358 (0) 9 5064140

Fax: 358 (0) 9 50641411

E-mail: finland@tainstruments.com

(continued on next page)

FRANCE

TA Instruments Division de Waters SA
1-3, Rue Jacques Monod
78280 Guyancourt
France
Phone: 33/1 30 48 94 60
Fax: 33/1 30 48 94 51
E-mail: france@tainstruments.com

GERMANY

Waters GmbH
UB TA Instruments
Helfmann – Park 10
D-65760 Eschborn
Germany
Phone: 49/6196 400600
Fax: 49/6196 4006060
E-mail: germany@tainstruments.com

INDIA

c/o Waters (India) Pvt. Ltd.
36A, II Phase
Peenya Industrial Area
Bangalore 560 058
India
Attn: M. Sundar
Phone: 91 80 28398963
Fax: 91 80 28398964
E-mail: India@tainstruments.com

ITALY

TA Instruments – Div. Di Waters S.p.A.
Via Achille Grandi, 27
20090 Vimodrone (Milano),
Italy
Phone: 39/02 2742 11
Fax: 39/02 250 1827
E-mail: italia@tainstruments.com

JAPAN

TA Instruments Japan
No. 5 Koike Bldg.
1-3-12 Kitashinagawa
Shinagawa-Ku, Tokyo 140-0001
Japan
Phone: 813 5479 8418 (Sales & Application)
Fax: 813 5479 7488 (Sales & Application)
Phone: 813 3450 0981 (For Service & Accounting)
E-mail: J-Marketing@tainstruments.com

MEXICO

TA Instruments – Waters LLC
Moras 822, Colonia Acacias
C.P. 03230 Mexico, D. F.
Phone: 52 55 52 00 18 60
Fax: 52 55 55 24 93 76
E-mail: mdominguez@tainstruments.com

THE NETHERLANDS

TA Instruments
A Division of Waters Chroma-tography bv
Postbus 379 / Florijnstraat 19
4870 AJ Etten-Leur
The Netherlands
Phone: 31/76 508 72 70
Fax: 31/76 508 72 80
E-mail: Netherlands@tainstruments.com
Shipping address:
Mon Plaisir 12
4870 AJ Etten-Leur
The Netherlands

NORWAY

Hvamstubben 17
N-2013 Skjetten
Norway
Phone: 47 63846055
Fax: 47 63846051
E-mail: norway@tainstruments.com

SPAIN

Waters Cromatografia, s.a.
Division TA Instruments
Parc Tecnologic del Valles
Ronda de Can Fatjo, 7-A
08290- Cerdanyola del Valles
Phone: 34 93 600 93 00
Fax: 34 93 600 93 60
E-mail: spain@tainstruments.com

SWEDEN

TA Instruments Division of Waters Sweden
Djupdalsvagen 12-14
PO Box 485
19124 Sollentuna
Claus Larsen
Phone: +46 (0) 8 555 11 521
Fax: +46 (0) 8 555 11 520
E-mail: sweden@tainstruments.com

TAIWAN

TA Instruments
Division of Waters Taiwan
Room 5, 4F., No23, Sec. 1
Chang-an E. Rd. Taipei 104 Taiwan
Phone: 886 (2) 2563 8880
Fax: 886 (2) 2563 8870
Email: TATaiwan@tainstruments.com

UNITED KINGDOM

TA Instruments Ltd
Block D, The Fleming Centre
Fleming Way, Royal Manor
Crawley, West Sussex, RH10 9NB
England
Phone: 44/1293 658900
Fax: 44/1293 658901
E-mail: uk@tainstruments.com or info@taeurope.co.uk