Q Series

Nitrogen Purge Cooler Accessory for DMA

Getting Started Guide

Notice

The material contained in this manual, and in the online help for the software used to support this instrument, is believed adequate for the intended use of the instrument. If the instrument or procedures are used for purposes other than those specified herein, confirmation of their suitability must be obtained from TA Instruments. Otherwise, TA Instruments does not guarantee any results and assumes no obligation or liability. TA Instruments also reserves the right to revise this document and to make changes without notice.

TA Instruments may have patents, patent applications, trademarks, copyrights, or other intellectual property covering subject matter in this document. Except as expressly provided in written license agreement from TA Instruments, the furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property.

TA Instruments Operating Software, as well as Module, Data Analysis, and Utility Software and their associated manuals and online help, are proprietary and copyrighted by TA Instruments. Purchasers are granted a license to use these software programs on the module and controller with which they were purchased. These programs may not be duplicated by the purchaser without the prior written consent of TA Instruments. Each licensed program shall remain the exclusive property of TA Instruments, and no rights or licenses are granted to the purchaser other than as specified above.

TA Instruments can accept no liability for loss or damage, however caused, arising from the faulty or incorrect use of its products. TA Instruments shall not be liable for any damages caused by interactions between exogenous materials (e.g. chemicals) and parts of the instrument. This includes interactions of gaseous, liquid or solid materials with parts of the instrument.

©2014 by
TA Instruments — Waters LLC
159 Lukens Drive
New Castle, DE 19720

Introduction

Important: TA Instruments Manual Supplement

Please click the [TA Manual Supplement](#) link to access the following important information supplemental to this Getting Started Guide:

- TA Instruments Trademarks
- TA Instruments Patents
- Other Trademarks
- TA Instruments End-User License Agreement
- TA Instruments Offices

Notes, Cautions, and Warnings

This manual uses NOTES, CAUTIONS, and WARNINGS to emphasize important and critical instructions. In the body of the manual these may be found in the shaded box on the outside of the page.

NOTE: A NOTE highlights important information about equipment or procedures.

CAUTION: A CAUTION emphasizes a procedure that may damage equipment or cause loss of data if not followed correctly.

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A WARNING indicates a procedure that may be hazardous to the operator or to the environment if not followed correctly.

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Safety

WARNING: The operator of this instrument is advised that if the equipment is used in a manner not specified in this manual, the protection provided by the equipment may be impaired.

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Instrument Symbols

The following labels are displayed on the Nitrogen Purge Cooler for your protection:

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Please heed the warning labels and take the necessary precautions when dealing with those parts of the instrument. The *DMA Nitrogen Purge Cooler Getting Started Guide* contains cautions and warnings that must be followed for your own safety.

Chemical Safety

WARNING: Do not use hydrogen or any other explosive gas in the NPC.

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WARNING: If you are using samples that may emit harmful gases, vent the gases by placing the instrument near an exhaust.

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WARNING: Air is not to be used as the purge gas with this cooler. Oxygen will be condensed into liquid and pushed into the DMA furnace.

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WARNING: MAY CAUSE SEVERE FROSTBITE. Can cause severe frostbite to the eyes and skin. Do not touch frosted pipes or fittings.

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WARNING: Do not overfill the Dewar. Leave an inch of space from the top of the Dewar when filling. Liquid nitrogen may boil over the sides when the Cooling Gas is turned on.

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Handling Liquid Nitrogen

This cooling accessory uses the cryogenic (low-temperature) agent, liquid nitrogen, for cooling. Because of its low temperature (-196°C [-321°F]), liquid nitrogen will burn the skin. When you work with liquid nitrogen, use the following precautions:

WARNING: Liquid nitrogen boils rapidly when exposed to room temperature. Be certain that areas where liquid nitrogen is used are well ventilated to prevent displacement of oxygen in the air.

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- 1 Wear a face shield, gloves large enough to be removed easily, and a rubber apron. For extra protection, wear high-topped, sturdy shoes, and leave your pant legs outside the tops.
- 2 Transfer the liquid slowly to prevent thermal shock to the equipment. Use containers that have satisfactory low-temperature properties. Ensure that closed containers have vents to relieve pressure.
- 3 The purity of liquid nitrogen decreases when exposed to air. If the liquid in a container has been open to the atmosphere for a prolonged period, analyze the remaining liquid before using it for any purpose where high oxygen content could be dangerous.

NOTE: Using the NPC poses all of the same risks as with any open Dewar of cryogenic liquid. Safe practices concerning the filling and transporting of small Dewars and use of Personal Protective Equipment is necessary.

The asphyxiant warning below applies to the use of liquid nitrogen. Oxygen depletion sensors are sometimes used where liquid nitrogen is in use.

**WARNING:
Potential Asphyxiant**

Liquid nitrogen can cause rapid suffocation without warning.

Store and use in an area with adequate ventilation.

Do not vent the Nitrogen Purge Cooler (NPC) in confined spaces.

Do not enter confined spaces where nitrogen gas may be present unless the area is well ventilated.

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Table of Contents

Introduction	3
Important: TA Instruments Manual Supplement	3
Notes, Cautions, and Warnings	4
Safety	5
Instrument Symbols	5
Chemical Safety	6
Table of Contents	9
Introducing the Nitrogen Purge Cooler	10
Overview.....	10
System Components	10
Accessory Specifications.....	11
Installing the Nitrogen Purge Cooler	13
Unpacking/Repacking the NPC.....	13
Preparing the Accessory	13
Inspecting the System	14
Choosing a Location	14
In	14
On	14
Near	14
Away from	14
Connecting the NPC to Your Q Series DMA.....	15
Setting up the NPC Accessory using Instrument Control Software.....	19
Filling the Dewar	20
Use and Maintenance	21
Using the NPC Accessory	21
Before You Begin	21
Running an NPC Experiment	22
Creating an Experimental Procedure	22
Editing an Experiment	23
Maintaining the Accessory	24
Cleaning the Accessory	24
Replacement Parts	24

Chapter 1:

Introducing the Nitrogen Purge Cooler

Overview

The Nitrogen Purge Cooler (NPC), is an optional accessory for extending the temperature range of the Q800 DMA standard furnace to -160°C. The NPC consists of a 2.5 L Dewar flask that contains a copper coil tube. The furnace is cooled by purging nitrogen gas through the copper coil immersed in the liquid nitrogen-filled Dewar flask. The NPC requires a nitrogen gas supply (25 to 120 psi) at a flow rate of 30 LPM, and access to a supply of liquid nitrogen and equipment necessary for safe handling, transportation, and pouring of the liquid.

NOTE: The NPC is for crash cooling and controlled heating only. Controlled cooling rates are possible using the Gas Cooling Accessory (GCA), or Air Chiller System (ACS-3).

System Components

The NPC Accessory has the following major hardware components:

- A Dewar flask with insulated lid, capable of holding approximately 2.5 L of liquid nitrogen, and with a precise length of copper tubing to adequately chill the nitrogen purge gas.
- A length of 1/8" OD tubing that connects the NPC to the Cooling Gas Outlet fitting on the rear of the DMA.

Accessory Specifications

The tables below detail the NPC Accessory technical specifications.

Table 1: NPC Accessory Characteristics

Dimensions of Dewar flask	Height: 15 in. Width: 12 in. Depth: 7 in. Volume: 2.5 L
Weight of Dewar flask	6 lbs (without liquid nitrogen)

Table 2: Accessory Cooling Gas Requirements

Gas	Nitrogen ONLY
Pressure	25–120 psi
Flow Rate	30 LPM

Table 3: Accessory Operating Environmental Conditions

Operating altitude	2000 meters maximum
Relative humidity	5% to 80% RH from 15°C to 31°C, decreasing to 66% RH at 35°C (non-condensing)

Table 4: Accessory Performance Specifications

Lowest temperature	-160°C
Time to reach low temperature	40 min
Time to reach -150°C	35 min
Heating rate: Minimum^a	3°C/min

a. With NPC off and dependent on ambient temperature.

Chapter 2:

Installing the Nitrogen Purge Cooler

Unpacking/Repacking the NPC

The instructions needed to unpack and repack the accessory are found as separate unpacking instructions in the shipping box and in the online documentation associated with the instrument control software. Be sure to read and perform the unpacking instructions prior to performing any procedures in this chapter.

Retain all of the shipping hardware and boxes from the accessory in the event you wish to repack and ship your accessory.

Preparing the Accessory

Before shipment, the NPC is inspected so that it is ready for operation upon proper installation. Only limited instructions are given in this manual; consult the online documentation for additional information. Installation involves the following procedures:

- Inspecting the accessory for shipping damage and missing parts
- Connecting the NPC tubes and fittings
- Filling the Dewar

CAUTION: To avoid mistakes, read this entire chapter before you begin installation.

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Inspecting the System

When you receive the NPC Accessory, look over the accessory and shipping container carefully for signs of shipping damage, and check the parts received against the enclosed shipping list.

- If the accessory is damaged, notify the carrier and TA Instruments immediately.
- If the accessory is intact but parts are missing, contact TA Instruments.

Choosing a Location

Choose a location for the accessory using the following guidelines. The NPC Accessory should be:

In

- A temperature-controlled area. Temperatures should be in the range of 20–35°C.
- A clean environment, preferably on the ground floor in the building.
- An area with ample working and ventilation space.

On

- A stable work surface.

Near

- Your Q Series DMA.

Away from

- Dusty environments.
- Exposure to direct sunlight.
- Direct air drafts (fans, room air ducts).
- Poorly ventilated areas.
- Noisy or mechanical vibrations.
- High traffic areas, where constant movements from passing personnel could create air currents or mechanical disturbances.

Connecting the NPC to Your Q Series DMA

To connect the NPC, access the rear panel of the DMA and follow the instructions below:

NOTE: Remove the foam insulation from the Dewar before making connections and filling the Dewar with nitrogen.

- 1 Place the Dewar next to the DMA.
- 2 Disconnect the GCA tube from the DMA.

Figure 1 Disconnect the GCA tube.

- 3 Connect the copper tubing to the GCA port on the DMA.

Figure 2 Connect copper tubing to the DMA.

- 4 Fit the tube insulation piece over the GCA port connection.

Figure 3 Tube insulation piece covering the GCA connection.

- 5 Connect nitrogen purge gas to the Cooling Gas Inlet port (see [Figure 4](#) below) on the DMA. Inlet pressure should be between 25 and 120 psi. Use nitrogen gas ONLY.
- 6 Remove the silicone tubing from the barb fitting on the Air Cool Outlet port, then remove the barb fitting.


Figure 4 Barb fitting with silicone tubing attached to the Air Cool Outlet port.

- 7 Place the Cooling Gas Outlet label directly over the Air Cool Outlet label so that “Air Cool Outlet” is completely obstructed.

- 8 Attach the quick-disconnect elbow to the Cooling Gas Outlet port.

Figure 5 Quick-disconnect elbow fitting attached to Cooling Gas Outlet port.

- 9 Insert the empty end of the 1/8" black tubing into the end of the elbow. Then screw the nut attached to the end of the tube to the Swagelok fitting on the copper tubing.

Figure 6 Black tubing (left) connected to the copper tubing in the Dewar (right).

NOTE: The 20" length of 1/8" tubing is critical. Do not shorten this tube. Shortening this tube can cause an excessive amount of nitrogen to be condensed into liquid and pushed into the furnace.

- 10 Proceed to the [next section](#) to set up the accessory in Thermal Advantage. Do not [fill the Dewar](#) until the next section has been completed.

Setting up the NPC Accessory using Instrument Control Software

CAUTION: Operation of this accessory is intended for the QNX-based DMA Q800 with instrument software version 21.2.73 or higher only.

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After successfully installing the NPC Accessory, access the DMA Thermal Advantage Software and follow the instructions below.

1 Select **Control** > **Select Accessory**


Figure 7 Select Accessory.

2 In the DMA Cooler Selection dialog box, select **NPC**.


Figure 8 Select NPC.

3 Proceed to [Filling the Dewar](#).

Filling the Dewar

Before the Dewar can be filled with liquid nitrogen, the copper tube must be purged with dry gas. This is necessary to prevent moisture from freezing in the copper tube when liquid nitrogen is added to the Dewar; the frozen moisture blocks the flow of the cooling gas.

- 1 Make sure the furnace is closed.
- 2 In Thermal Advantage, select **Cooling Gas > On** from the **Control** tab.


Figure 9 Turn on Cooling Gas in Thermal Advantage software.

- 3 Allow the dry gas to purge through the copper tube for at least 20 minutes, or until the copper tube is dry.
- 4 Turn the Cooling Gas off (the same way it was turned on in the software).
- 5 Fill the Dewar with liquid nitrogen.

WARNING: Do not overfill the Dewar. Liquid nitrogen may boil over the sides when the Cooling Gas is turned on.

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Use and Maintenance

Using the NPC Accessory

All of your NPC Accessory experiments will follow the same general outline. In some cases, not all of these steps will be performed. The majority of these steps are performed using the instrument control software. See the instrument control online help for instructions on performing these actions.

- Filling the Dewar
- Selecting and preparing the sample
- Creating or choosing a test procedure and entering experiment information through the TA Instruments instrument control software
- Loading the prepared sample
- Starting the experiment

Please note the following when performing an experiment:

- One full Dewar (approximately 2.5 L of liquid nitrogen) will typically be consumed for each experiment.
- A full Dewar is sufficient for performing the following method when the starting temperature is 40°C:
 - External Event: ON
 - Equilibrate: -150°C
 - Isothermal: 5 minutes
 - External Event: OFF
- When the NPC is connected to a bottle of dry, compressed nitrogen at approximately 2000 psi, seven complete experiments as shown above can be performed.

Before You Begin

Before using the NPC Accessory, ensure that the DMA Q800 Dynamic Mechanical Analyzer is installed properly. Also make sure you have:

- Connected the NPC
- Powered on the DMA
- Filled the Dewar as necessary
- Specified the NPC in the instrument control software
- Become familiar with controller operations

It is a good idea to purge the copper coil with dry nitrogen just prior to filling the Dewar with liquid nitrogen. Select **Control > Cooling Gas > On** to purge the coil with dry nitrogen. After a few minutes, select **Control > Cooling Gas > Off** to turn the purge gas off again.

Running an NPC Experiment

Creating an Experimental Procedure

- 1 From the **Procedure** tab, select the test you want to run from the **Test** drop-down menu

Summary Procedure Notes

Procedure Information

Test: Temp Ramp / Freq Sweep

Notes: Material is heated at a constant rate. While heating, the material is deformed (oscillated) at a constant amplitude (strain) over a range of frequencies and the mechanical properties measured.

Temperature Ramp / Single Frequency

Amplitude: 15.0000 μm

Strain: 0.0000 %

Advanced...

Post Test...

Start temperature: -150.00 °C

Soak time: 5.00 min

Final temperature: 150.00 °C

Ramp rate: 3.00 °C/min

Hold time at final temperature: 30.00 min

Method / Frequency Table

Figure 10 Procedure tab > Test menu.

- 2 Change the ramp/frequency and temperature fields as needed.
- 3 The experiment steps display on the right in the **Running Segment Description** section. Click **Apply** when you have the experiment set up to your satisfaction.

#	Running Segment Description
1	Data storage Off
2	Equilibrate at -150.00 °C
3	Isothermal for 5.00 min
4	Data storage On
5	Ramp 3.00 °C/min to 150.00 °C

Figure 11 Experiment steps.

Editing an Experiment

- 1 In the **Procedure** tab, select **Custom** from the **Test** drop-down menu, and then click **Editor**.


Figure 12 Custom > Editor.

- 2 A **Method** dialog box displays. Drag experiment steps from the Segment list and drop them into the Segment Description section where needed. Click **OK** when finished.


Figure 13 Edit the experiment steps.

NOTE: The segment “External Event On/Off” refers to turning the Nitrogen Purge Cooler on or off.

Maintaining the Accessory

The primary maintenance procedures described in this section are the customer's responsibility. Any further maintenance should be performed by a representative of TA Instruments or other qualified service personnel.

Cleaning the Accessory

To clean the NPC Accessory, wipe down the exterior of the Dewar with a damp, soft cloth.

CAUTION: Do not use harsh chemicals, abrasive cleansers, steel wool, or any rough materials to clean the cabinet, as you may scratch the surface and degrade its properties.

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Replacement Parts

Table 5: Replacement Parts for the NPC Accessory

Part Number	Description
986316.901	Insulation Replacement Kit NPC